

INSTRUCTIONS FOR USE

Transform S MicroPath

Nickel-titanium glide path and conservative preparation files

FOR PROFESSIONAL DENTAL USE**01 Intended use and product range**

For enlargement of an established glide path and limited canal preparation with the designated file variants, by dental professionals trained in endodontics. Identify the file by its article number, tip size, taper and length.

Glide path files	13/.03, 15/.03 and 17/.03; 21, 25, 29 mm
Short variant	15/.05, 17 mm
Extension files	20/.03 and 25/.03; 21, 25, 29 mm

Not all size–length combinations are available. Range and final device specifications require approval.

02 Composition and presentation

Nickel-titanium working part with a handpiece shank and depth stop. This draft covers the planned sterile, single-use presentation. Final component materials and coatings require confirmation.

03 Contraindications and safety

Do not use in patients allergic to a constituent material. Nickel-titanium contains nickel. Avoid powered instrumentation in severe, abrupt apical curves where safe progression cannot be achieved.

Use rubber dam isolation, eye protection and standard infection-control precautions. Sharp instruments can cause injury. Misuse or damage may cause fracture, canal transportation, perforation or cross-contamination.

Do not use expired products, open or damaged sterile packs, or files with corrosion, cracks, unwound flutes or other damage. Do not reuse or resterilise single-use files.

Use minimal pressure. Stop when a file binds or ceases to advance. Inspect throughout treatment; fatigue may occur without visible damage.

04 Preparation and operating mode

Continuous clockwise rotation with a suitable torque-controlled endodontic motor.

Before issue: confirm speed, torque and the program for each file group. Do not use the draft protocol clinically.

Assess clinical findings, radiographs and canal anatomy. Prepare access and determine working length. Open the pack aseptically before use; check file identity and stop position. Confirm secure handpiece retention and the approved program.

05 Clinical sequence

1. Locate and explore the canal with small hand files. Establish working length and a reproducible path before introducing any powered glide path instrument.
2. For the base sequence, progress through 13/.03, 15/.03 and 17/.03 only as permitted by the approved protocol and canal anatomy. Do not skip to a larger file when the smaller instrument cannot advance passively.
3. Advance using light, short movements. On resistance, withdraw and clean the file, irrigate and recapitulate manually. Reconfirm the path before proceeding.
4. The 15/.05, 17 mm short variant has a different geometry. Its validated depth and role must be specified before use; it is not automatically interchangeable with the .03 glide path instruments.
5. The 20/.03 and 25/.03 extensions are a separate preparation option. Use only after the required path has been achieved and their dedicated settings and sequence have been approved. Reassess apical size before further shaping.

06 Irrigation and instrument care

Maintain irrigation and remove debris during preparation. Clear flutes without damaging edges. Do not immerse the handle or entire instrument in irrigant. Before issue, confirm compatible solutions, concentrations and contact times for this device.

07 Storage disposal and reporting

Protect unopened packs from moisture, sunlight and damage. Before issue, confirm temperature, humidity and shelf-life limits. Follow the final pack expiry date and batch information. Discard used files in an approved sharps container under the practice's clinical-waste procedure.

Record article and batch details where practicable. Report defects or adverse incidents to the NZ distributor, retaining packaging and relevant device details.

08 New Zealand distributor

ENDOTECHNZ LIMITED

3c/14 Beaumont St, Auckland 1010, New Zealand
Steve Shepherd · 021 543 414