

INSTRUCTIONS FOR USE

Transform S RG

Nickel-titanium reciprocating root canal preparation system

FOR PROFESSIONAL DENTAL USE**01 Intended use and product range**

For glide path enlargement and mechanical shaping of root canals with the designated instruments, by dental professionals trained in endodontics. Identify the file by its article number, tip size, taper and length.

Shaping files	Small 20/.07; Primary 25/.07; Medium 35/.06; Large 45/.05
Glider	15/.02; 21, 25 and 31 mm
Shaping lengths	21, 25 and 31 mm
Orifice shaper	18/.04, 16 mm

Not all size-length combinations are available. Range and final device specifications require approval.

02 Composition and presentation

Nickel-titanium working part with a handpiece shank and depth stop. This draft covers the planned sterile, single-use presentation. Final component materials and coatings require confirmation.

03 Contraindications and safety

Do not use in patients allergic to a constituent material. Nickel-titanium contains nickel. Avoid powered instrumentation in severe, abrupt apical curves where safe progression cannot be achieved.

Use rubber dam isolation, eye protection and standard infection-control precautions. Sharp instruments can cause injury. Misuse or damage may cause fracture, canal transportation, perforation or cross-contamination.

Do not use expired products, open or damaged sterile packs, or files with corrosion, cracks, unwound flutes or other damage. Do not reuse or resterilise single-use files.

Use minimal pressure. Stop when a file binds or ceases to advance. Inspect throughout treatment; fatigue may occur without visible damage.

04 Preparation and operating mode

Reciprocating shaping. Use a validated motor program and motion parameters.

Before issue: confirm speed, torque and the program for each file group. Do not use the draft protocol clinically.

Assess clinical findings, radiographs and canal anatomy. Prepare access and determine working length. Open the pack aseptically before use; check file identity and stop position. Confirm secure handpiece retention and the approved program.

05 Clinical sequence

1. Explore the canal manually and establish working length and a reproducible path before mechanical preparation. Enlarge coronal access with the designated orifice shaper only where indicated.
2. Use the glider to enlarge the established path according to its approved motion and settings. Do not assume the glider, shaper and orifice instrument use the same program.
3. Select the shaping instrument according to the approved RG selection protocol and canal anatomy. The Primary 25/.07 is the reference starting option under review; use of Small 20/.07 when progression is restricted requires a patent, reproducible path.
4. Use light, short advancing movements. Remove the file frequently for cleaning and inspection; irrigate and recapitulate. Stop if progression ceases and reassess before re-entry.
5. Gauge the apical preparation. Medium 35/.06 or Large 45/.05 may be required for wider canals under the approved protocol; do not enlarge automatically through every size.

06 Irrigation and instrument care

Maintain irrigation and remove debris during preparation. Clear flutes without damaging edges. Do not immerse the handle or entire instrument in irrigant. Before issue, confirm compatible solutions, concentrations and contact times for this device.

07 Storage disposal and reporting

Protect unopened packs from moisture, sunlight and damage. Before issue, confirm temperature, humidity and shelf-life limits. Follow the final pack expiry date and batch information. Discard used files in an approved sharps container under the practice's clinical-waste procedure.

Record article and batch details where practicable. Report defects or adverse incidents to the NZ distributor, retaining packaging and relevant device details.

08 New Zealand distributor

ENDOTECHNZ LIMITED

3c/14 Beaumont St, Auckland 1010, New Zealand
Steve Shepherd · 021 543 414