

INSTRUCTIONS FOR USE

Transform S C Pilot

Stainless steel C Pilot hand exploration files

FOR PROFESSIONAL DENTAL USE**01 Intended use and product range**

For manual canal scouting and initial glide path development in narrow canals, by dental professionals trained in endodontics. Identify the file by its article number, tip size, taper and length.

Tip sizes	6, 8, 10, 12.5 and 15
Profile	C Pilot profile; final taper specification to confirm
Lengths	19, 21 and 25 mm

Not all size–length combinations are available. Range and final device specifications require approval.

02 Composition and presentation

Stainless steel working part with a manual handle and depth stop. This draft covers the planned sterile, single-use presentation. Final component materials and coatings require confirmation.

03 Contraindications and safety

Do not use in patients allergic to a constituent material. Check the alloy declaration for patients with metal allergies. Do not force a file through a blocked or unidentifiable path.

Use rubber dam isolation, eye protection and standard infection-control precautions. Sharp instruments can cause injury. Misuse or damage may cause fracture, canal transportation, perforation or cross-contamination.

Do not use expired products, open or damaged sterile packs, or files with corrosion, cracks, unwound flutes or other damage. Do not reuse or resterilise single-use files.

Use minimal pressure. Stop when a file binds or ceases to advance. Inspect throughout treatment; fatigue may occur without visible damage.

04 Preparation and operating mode

Manual use only. Do not connect to a motor.

Before issue: confirm manual technique and movement limits.

Assess clinical findings, radiographs and canal anatomy. Prepare access and determine working length. Open the pack aseptically before use; check file identity and stop position. Maintain a controlled manual grip.

05 Clinical sequence

1. Assess canal anatomy and choose a small file appropriate to the canal entrance. Set the stop to the planned exploration depth.
2. Explore manually using light, controlled movements. Withdraw frequently to remove debris, inspect the instrument and irrigate.
3. Do not force advancement into an obstruction. Reassess the canal direction and access when the file does not progress.
4. Increase size gradually as a reproducible path develops. Confirm working length before full-length glide path preparation.
5. Before changing to powered files, verify that the required glide path for that system has been achieved. Discard instruments showing damage or loss of cutting integrity.

06 Irrigation and instrument care

Maintain irrigation and remove debris during preparation. Clear flutes without damaging edges. Do not immerse the handle or entire instrument in irrigant. Before issue, confirm compatible solutions, concentrations and contact times for this device.

07 Storage disposal and reporting

Protect unopened packs from moisture, sunlight and damage. Before issue, confirm temperature, humidity and shelf-life limits. Follow the final pack expiry date and batch information. Discard used files in an approved sharps container under the practice's clinical-waste procedure.

Record article and batch details where practicable. Report defects or adverse incidents to the NZ distributor, retaining packaging and relevant device details.

08 New Zealand distributor

ENDOTECHNZ LIMITED

3c/14 Beaumont St, Auckland 1010, New Zealand
Steve Shepherd · 021 543 414