

Transform S™ ET

Compatible with X7 and HyFlex

Establish the glide path, then select .04 or .06 taper.
ENDOTECH ROTARY · 400-500 rpm · 2 N·cm · MicroPath 13/.03: 1 N·cm

ENTRY CHECK

Assess anatomy, isolate with a rubber dam, scout with a small hand file and establish working length (WL).

THE INSTRUMENTS



ORIFICE
SHAPER
17/.08 · 17 mm

Available in sizes 15 to 45.

SEQUENCE + CLINICAL DECISIONS

- 1 Glide path**
MicroPath 15 or 17; 13 if needed
After hand scouting, select 15/.03 or 17/.03 for routine glide preparation; use 13/.03 for complex or narrow paths. Confirm a smooth, reproducible path to WL.
- 2 Open the entrance**
17/.08 · 17 mm
Use the dedicated rotary opener only in the established coronal path. Brush away from the furcation; stop before curvature and resistance.
- 3 Initial preparation**
17/.04
Follow the glide path with light pressure, removing frequently to clean and inspect. Reconfirm WL after coronal enlargement.
- 4 Develop the shape**
25/.04; 20/.04 if needed
Try 25/.04 passively to WL. If it cannot advance, use 20/.04, then return to 25/.04. Never drive through resistance.
- 5 Choose the endpoint**
Gauge size and taper
Select further enlargement from the apical assessment and anatomy. Use 25/.06 or larger .04/.06 files only when justified.
- 6 Complete treatment**
Clean → cone fit → seal
Record final size/taper and WL. Complete disinfection, check the selected cone at WL, then obturate under length control.

GLIDE → OPEN → 17/.04 → 25/.04 → GAUGE

Optional 20/.06 preflare

EndoTech suggestion: 20/.06 may preflare toward mid-root only within a confirmed, sufficiently wide straight coronal path. Stop short of curvature and resistance; mid-root is not a mandatory depth.

Inspect between instruments

Withdraw. Clean and inspect flutes. Irrigate, recapitulate with a small hand file and re-irrigate. If the path is no longer reproducible, rebuild it before advancing.

Advantage of .06 preflare

Taking a size 20/.06 file to mid-root before 17/.04 and 25/.04 can improve cone fit in the final preparation. Verify the final cone at working length.

Transform S™ ET

Compatible with X7 and HyFlex

Choose the route that fits the canal.

01

ET compatibility

Compatible with X7 and HyFlex. This card uses the X7-style 17/.04 to 25/.04 technique with MicroPath glider, opener and optional preflare stages. Use EndoTech motor settings as printed below.

03

Orifice opener

The orifice opener is a unique 17/.08, 17 mm file. It was created specifically for endodontists and dentists who find larger-tipped orifice shapers too large to easily enter smaller calcified orifices, and gives far better positioning for redefining the orifice in more difficult canals.

02

Alternative crown down

ET files can equally be used successfully in a crown down technique: 25 → 20 → 17, in .04 or .06 taper.

04

Instructions for use

Read the instrument designation on the packaging before use. ET shaping files are available in sizes 15 to 45. The orifice shaper is 17/.08, nominal 17 mm. The pictures are not a size or length scale.

MOTOR SETUP

ROTARY: 400-500 rpm. ET, the dedicated 17/.08, 17 mm opener and MicroPath 15/.03 / 17/.03: 2 N-cm. MicroPath 13/.03: 1 N-cm.

Orifice opener · 17/.08 · 17 mm



Transform S™ ET

MicroPath glide followed by ET opening and shaping.
ENDOTECH ROTARY · 400-500 rpm · 2 N·cm

Compatible with X7 and HyFlex

15/.03 → 20/.06 → 17/.04 → 25/.04



Use 20/.06 towards mid-root only within the established, sufficiently wide straight coronal path. Stop before curvature or resistance. Clean, inspect, irrigate and recapitulate between instruments.