



SHAPING SYSTEMS

TransformSTM

PT + RG

PROGRESSIVE TAPER ROTARY
COMPARE WITH PROTAPER GOLD

RECIPROCATING SHAPING
COMPARE WITH WAVEONE GOLD

Two shaping motions.
One TransformS family.

Choose the sequence that defines your
shaping workflow.



PT - F2



RG - PRIMARY

● PT progressive taper
Compare with ProTaper Gold

● RG reciprocating
Compare with WaveOne Gold

TransformTM PT

Progressive taper.

Shape, then finish with an exceptionally versatile file.

01. CORONAL



Dedicated orifice opener.

SX · 19 mm

02. SHAPING



Progressive taper shaping
(S1-S2) files available in
three working lengths.

S1 · S2
21 / 25 / 31 mm

03. FINISHING



Finishing files (F1-F5)
available in three working
lengths.

F1 · F2 · F3
21 / 25 / 31 mm

COMPARE WITH PROTAPER GOLD
AND EDGETAPER PLATINUM

SINGLE SIZES AND ASSORTED
SIZES AVAILABLE IN
21 / 25 / 31 MM LENGTHS

FOR PROFESSIONAL DENTAL USE.

TransformTM RG

Reciprocating shaping.

Opening, glide path and four shaping sizes.

01. OPENING



Dedicated orifice shaper.

ORIFICE OPENER

16 mm · Pack of 3

02. GLIDE PATH



Glide-path support for the
RG shaping sequence.

GLIDE PATH

21 / 25 / 31 mm · Packs of 3 or 6

03. SHAPING



Four shaping sizes, with
individual and assorted pack
options.

SMALL · PRIMARY

MEDIUM · LARGE

21 / 25 / 31 mm

COMPARE WITH
WAVEONE GOLD

PACKS OF 3 OR 6
ASSORTED PACK OF 4

FOR PROFESSIONAL DENTAL USE.

SYSTEM SELECTION

Two shaping pathways.

Select by motion, sequence and case requirements.

PT PROGRESSIVE TAPER Continuous rotary		RG RECIPROCATING Reciprocating shaping	
Sequence	SX → S1 → S2 → finishing file	Sequence	Orifice shaping → glide path → shaping
Finishing range	F1 through F5	Shaping range	Small · Primary · Medium · Large
Lengths	19 mm SX · 21 / 25 / 31 mm range	Lengths	16 mm opener · 21 / 25 / 31 mm range
Individual pack	6 files	Individual pack	3 or 6 files



TransformS™
PT + RG

Two shaping motions.
One clear family.

Incredible FlexION technology at incredible value.